

## In Situ Cell Death Detection Kit, Fluorescein

Cat. No. 11 684 795 001 50 tests

## In Situ Cell Death Detection Kit, TMR red


Cat. No. 12 156 792 001 50 tests

<b>Type</b>	Direct TUNEL labeling assay
<b>Useful for</b>	Detection of DNA strand breaks in apoptotic cells by flow cytometry or fluorescence microscopy
<b>Samples</b>	Cells in suspension, adherent cells, cell smears, frozen or paraffin-embedded tissue sections
<b>Method</b>	End-labeling of DNA with fluorescein-dUTP or tetramethylrhodamine-dUTP (TMR-dUTP), followed by direct analysis of fluorescent cells
<b>Time</b>	1–2 h (+ sample preparation, permeabilization, etc.)

**Significance of kit:** This two *In Situ* Cell Death Detection Kits, measure and quantitate cell death (apoptosis) by labeling and detection of DNA strand breaks in individual cells by flow cytometry or fluorescence microscopy. The kits offer a direct TUNEL detection method, for maximum sensitivity and minimal background.

**Test principle:** The assays use an optimized terminal transferase (TdT) to label free 3'OH ends in genomic DNA with fluorescein-dUTP or TMR-dUTP. The procedure involves:

- 1 Fixing and permeabilizing apoptotic cells.
- 2 Incubating the cells with the TUNEL reaction mixture containing TdT and fluorescein-dUTP or TMR-dUTP. During this incubation step, TdT catalyzes the attachment of fluorescein-dUTP or TMR-dUTP to free 3'OH ends in the DNA.
- 3 Visualizing the incorporated fluorescein with a flow cytometer and/or a fluorescence microscope (fluorescein/TMR red).




▲ Figure 24: Schematic showing the principle of the *In Situ* Cell Death Detection Kits, Fluorescein and TMR red.

For a detailed overview of the steps in the procedure, see Flow Chart 7.

**Sensitivity:** The enzymatic labeling allows the detection of an apoptotic event that occurs, prior to changes in morphology and even before DNA fragments become detectable in the cytoplasm<sup>22</sup>. It detects early stage of DNA fragmentation in apoptotic cells. This is especially important if apoptosis is studied *in vivo*, e.g., in tissue sections, since apoptotic cells are rapidly and efficiently removed *in vivo*.

**Specificity:** The amount of DNA strand breaks in apoptotic cells is so large that the degree of cell labeling in these assays is an adequate discriminator between apoptotic and necrotic cells<sup>19</sup>.



▲ Flow Chart 7: Assay procedure, *In Situ* Cell Death Detection Kits (Fluorescein or TMR red).

**Can be used to assay:**

- Cells in suspension (permanent cell lines, normal and tumor cells *ex vivo*)
- Cytopspins, cell smears
- Adherent cells cultured on chamber slides
- Frozen tissue sections
- Formalin-fixed, paraffin-embedded tissue sections

**Kit contents**


1. Enzyme solution (TdT), 5 tubes
2. Labeling solution (nucleotide mix), 5 tubes

**Typical results:** See Figures 25–28.


**Technical tips:** For more information on the use of the kit for flow cytometric analysis, see page 119 in the Appendix of this guide.

**Other applications:** For more examples of how the *In Situ* Cell Death Detection Kits (Fluorescein or TMR red) can be used in the lab, see Appendix, pages 134–136.

# A




▲ **Figure 25: Detection of apoptotic cells by flow cytometry using the *In Situ* Cell Death Detection Kit, Fluorescein.** HL60 cells were cultured in the absence (A) or presence (B) of 2 µg/ml Camptothecin for 3 h at 37°C. Cells were incubated either with TUNEL reaction mixture (■) or label solution as negative control (■) or PBS for autofluorescence (□).




◀ **Figure 26: Detection of apoptotic cells (green) by fluorescence microscopy in a tissue section from rat.** A tissue section from a rat spinal cord was prepared and assayed with the *In Situ* Cell Death Detection Kit, Fluorescein. The treated section was viewed under a fluorescence microscope. (Photomicrograph was kindly provided by R. Gold, University of Würzburg, Germany.)

**Result:** A subpopulation of apoptotic cells, scattered throughout the tissue section, are intensely stained (green) by the TUNEL treatment and are easily visible under the microscope.



◀ **Figure 27: Cell suspension stained with the *In Situ* Cell Death Detection Kit, Fluorescein.** U937 cells induced with 4 µg/ml camptothecin, showing positive staining of apoptotic nuclei.

**Note:** This figure shows a high number (>80%) of apoptotic cells. To avoid detecting cells that are undergoing secondary necrosis, analyze cells earlier in the process after induction of apoptosis.



◀ **Figure 28: Rabbit endometrium, stained with the *In Situ* Cell Death Detection Kit, TMR red, and viewed under a fluorescence microscope.** Apoptotic nuclei stain bright red, limited fluorescence is visible in background tissue.